




30MA REAL PORTABLE X-RAY MACHINE

30mA-85KVP Full
Wave Rectified X-Ray
Generator having
Stationary Anode X-Ray
Tube as per ISI
Specifications
Truly Portable Carry
away box type machine.



Technical Details:



Features



Voltage Overload Protection.



Exposure Time can be set by Soft Button.



Hand Switch for Exposure in single touch.



Machine Safety Devices - Electronic Sensor to protect X-Ray Tube & High Voltage transformer from overloading, Quick Trip Overload Circuit Breaker, Fuses for all circuits in the machine



Control

(with inbuilt Tube Head) - Compact control panel with soft touch switch having independent Voltmeter to indicate line voltage. mA meter, Booster transformer is provided in control panel for stabilization of filament voltage. Technique selector for selection of Radiographic mA, kvp selector from 30 to 85 & Time Selector for selection of Radiation time, Standby & Exposure switches.


Accessories

Power Supply Cable,
Exposure switch with cable.
LCD collimator





Contact Us



SPOX HEALTHCARE PVT. LTD.

B-15 Sector-6, Noida, Uttar Pradesh 201301 (INDIA)

+91 98113 14143 www.spox.in

Certified & Approved by:

