



PRESENTING NEW RANGE OF X-RAY RADIATION PROTECTION ACCESSORIES



Lightweight and protective lead and lead-free radiation protection aprons for all types of procedures. Choose an apron and customize with color, size and optional accessories.

Some More Color and Designs



SX-LA-G



SX-LA-R



SX-LA-RB



SX-LA-P




SX-LA-BY

X Ray Radiation Protective Lead Apron

Wearing this radiation protection lead aprons greatly reduces an individual's exposure to radiation. The regular style contains 100% lead making them the heaviest apron. Lead composite aprons use a mixture of lead & other lighter weight metals making these lead aprons are lighter than regular lead aprons.

COLOURS

- 1 Sea Green
- 2 Parrot Green
- 3 Pink
- 4 Red
- 5 Yellow
- 6 White




Lead Apron SX-LANS-0.25mm Pb Lead Equivalent
Lead Apron SX-LANS-0.5mm Pb Lead Equivalent


Specifications:-

- Every Revolution Apron comes with Cool Wear moisture-wicking fabric
- Partial overlap to increase protection
- Two-tone color blocking
- In revolution apron we have cool wear moisture-wicking fabric
- Professionally designed for an exceptional fit and comfort
- Each color has been custom created for the revolution line and it is lightweight, waterproof and abrasion resistant.
- Designed for ease of wearing
- Adjustable back strap for snug fit.


SOME MORE COLOR AND DESIGNS




SX-LANS-W




SX-LANS-Y



SX-LANS-SG



SX-LANS-R



SX-LANS-G



SX-LANS-O



SX-LANS-P

CE & FDA Approved Lead Apron



CE & FDA Lead Apron 0.35mm Pb Lead Equivalent



CE & FDA Lead Apron 0.5mm Pb Lead Equivalent

X-Ray Protection Thyroid Collar

This thyroid collar provide protection for the neck and thyroid area from radiation exposure, its velcro-belt offer the flexibility to fit any neck size. It helps to cut-back patient exposure throughout intraoral radiography. Light-weight and comfortable.



X-Ray Protection Lead Gloves

Comfortable fit with maximum dexterity. Soft removable, washable and replaceable inner liners* reduce odors and permit frequent laundering. Expanded durable vinyl outer covers are removable for inspection of lead core and permit replacement when they become worn. The lead vinyl cores are precision molded for consistent 0.5 mm lead equivalency and long life.



Features:

- 0.25mm Pb
- Tough, durable vinyl outer coverings
- Soft, Flexible and Comfortable
- Removable, washable and replaceable Cool Wear moisture-wicking inner liners

More Colour



SX-LG-01



SX-LG-02



SX-LG-03



SX-LG-04



SX-LG-05



SX-LG-06

X-ray Radiation Protective Lead Cap

Specification

- One Size Fits All
- Exclusive Revolution outer fabric
- Scatter radiation protection
- Cool Wear Moisture-Wicking Fabric
- Hook & loop back closure for customizable fit




X-ray Radiation Protective Lead Goggles

Specification

- Light and comfortable molded frames for easy wearing. It covers front and side.
- Lead equivalency : 0.5mm
- Specially designed to suit over the eyes, it is a wonderful selection for those requiring smart splash protection from radiation.
- The large, clear facet shields defend the user from lateral exposure whereas eliminating the "tunnel" effects of look-alikes.



X-ray Radiation Protective Lead Mask



This radiation protecting mask eliminates blind spots caused by normal glasses and provides you a large space of protection from radiation. Great for people who wear glasses or who need to shield over simply their eyes. The leaded acrylic has a 0.1mm lead equivalency. It provides optimum performance for those people seeking a lead face protect providing a bigger space of protection than lead glasses. These liberally sized radiotherapy masks give maximum-width coverage to forestall exposure to scatter radiation from all angles. The adjustable x-ray face shields are often positioned at any angle to dispel scattered radiation waves.

SPOX HEALTHCARE PVT. LTD.

B-15 Sector-6, Noida, Uttar Pradesh 201301 (INDIA)

+91 98113 14143 | www.spox.in

Certified & Approved by:

