



Excellence Comes Standard



Driving The Digital Radiology Industry




ABOUT

SPOX Healthcare is a leading manufacturer of medical diagnostic X-Ray machines, specializing in producing Radiology and Imaging product for Human and Veterinary Healthcare. Combining vast experience, technical expertise and sound business strategies, we have steadily grown over the years and have demonstrated potential to become a major provider for Diagnostic Imaging Products, its spares, accessories and Radiation Protection Equipments. All equipments are vigorously tested to maintain highest standards of quality. SPOX has an ever expanding product range to meet changing customer needs worldwide. From its position of a leader in India, SPOX is spreading its wings wider, reaching out to distant and demanding markets all over the world. Our X-Ray machines are working to the satisfaction of our users who are in Government and Private Sectors.

PRODUCT RANGE


- Handheld Portable X-RAY
- Integrated Digital Mobile X-Ray
- DR system
- X-Ray Machine 30kW 40 kW 50kW
- C-ARM
- Radiation Protection Accessories
- Dental X-Ray Machine
- Film & Printer
- Spare Parts



Digital X-Ray with AI



WE are delighted to announce that a renowned film and television actor, **Sanjay Swaraj**, has been appointed as the **brand ambassador** of SPOX HealthCare Private Limited in India. This collaboration marks a significant milestone in our journey, reflecting our commitment to excellence and innovation in the healthcare sector.



HANDHELD PORTABLE X-RAY

4KG X-RAY MACHINE



GENERATOR TYPE	Line Frequency
mAs RANGE	.01-30sec
kV RANGE	70 kVp
mA RANGE	10 mA
FOCAL SPOT	0.8 mm
ANODE TYPE	stationary
POWER INPUT	190-230 V
DISPLAY	LCD
BODY TYPE	FIBER
STAND (Optional)	TRIPOD or SCISSOR ARM

30mA Real Portable X-RAY MACHINE



- 1.8mm - Focal Spot
- Weight - 15 kg
- 30mA - 85kVp

50mA Real Portable X-RAY MACHINE



- 1.8mm - Focal Spot
- Weight - 22 kg
- 50mA - 85 kVp

PORTABLE (HANDHELD) X-RAY

MODEL: SPOX SMART-5HS



Upto 80kV

Improved image penetration
> Improved bone identification

Upto 5mA

Reduced required irradiation time
> Reduced motion blurring

Low Dose

Minimal radiation exposure

High Resolution

Provides detailed and clear imaging

0.8mm

Focal Spot
> High resolution

400W

Power Output
> Ensuring Precise Results

Lightweight 3.5 kg

Easy to carry and transport

Wireless Remote Switch (Optional)

Offers remote control convenience.

INTEGRATED DIGITAL MOBILE X-RAY

4kW 6kW 8kW



MODEL	SPOX Ultimate 4kW	SPOX Ultimate 6kW
POWER OUTPUT	4KW	6KW

X-RAY GENERATOR & TUBE

TYPE	HF GENERATOR
mAs RANGE	0.1-320mAs
kV RANGE	40-110kV
mA RANGE	5-100mA/150mA;150mA@40kV/150mA/100ms
FOCAL SPOT	1.8mm
ANODE TYPE	stationary

SPOX IGZO DETECTOR



MODEL	SPOX 43i	SPOX 35i
APPLICATION	General Radiography	
CONVERSION SCREEN(SCINTILLATOR)	CsI/GOS	
TFT PANEL	a-Si/IGZO	
DIMENSIONS	461 x 461 x 15mm	385 x 461 x 15mm
WEIGHT INCL. BATTERY	<3.9Kg	<3.3Kg
PIXEL PITCH	140µm	
IMAGE DATA (AD CONVERSION)	16bit	
WIRED/ WIRELESS COMMUNICATION	Gigabit Ethernet/IEEE 802.11 ac	
TRIGGER MODE	High Sensitive 2XAED®	
DURABILITY	180Kg upto 50cm IP55	
*WEIGHT BEARING FOR WHOLE SURFACE *DROP TEST WATER & DUST RESISTANCE		
INNER STORAGE	<200 Images	<150 Images
BATTERY TYPE	Lithium-Ion-Polymer	
BATTERY'S FULL OPERATION TIME(MINIMUM)	Battery 12 hours	Battery 14 hours
IMAGE AREA SIZE	3,072 x 3,072mm	2,500x 3,052mm
X-RAY VOLTAGE MODE	40~150 kVp	
CALIBRATION MODE	ACC(Auto Calibration Control)/ Manual	
VOLTAGE	AC100~240V, 50/60 Hz, DC18V 3.5A	
ENVIRONMENT	Storage:- 10°C~60°C, 10% ~ 90% RH (non-condensing) Operation: 10°C~35°C, 20% ~ 75% RH (non-condensing)	

HF X-Ray MACHINE

30kW- 40kW- 50kW



C-ARM



LIVE IMAGE VIEWING & COMMUNICATION

- 10" touch console
- Real-time image support
- Easy workflow and operations



High power output
5 kw



Touch Monitor
Main: 27" or 2X19" (optional)
Sub: 10" console



Detector 9" X 9"
High resolution FPD



DICOM Services
DICOM 3.0
Connectivity to PACS

Radiation Protection Accessories



**Light Weight Vinyl Lead Apron X-Ray Vest Protection
&
Lead Free Apron (Ultra Light Weight Vinyl)**



**X-Ray Radiation Protective
Lead Goggles
0.5mm (Tested In NABL Lab)**



**X-Ray Radiation Protection
Lead Face
Shield 0.12mm Lead Pb
(Tested In NABL Lab)**



**X-Ray Radiation Protection
Lead Cap
0.5mm/0.25mm
(Tested In NABL Lab)**



**X-Ray Radiation Protective
Lead Goggles
0.75mm (Tested In NABL Lab)**



**X-Ray Protection Lead Thyroid Collar
Lead Equivalent (Tested In NABL Lab)**



**X-Ray Radiation Protection
Lead Gloves**



Lead Gonad Shield

Dental X-Ray Machine



Wall Model



**Scissor Arm
Floor Model**



XRDC70



**Scissor Arm
Wall Model**



Floor Model

PRINTERS



**LASER
PRINTER**



**THERMAL
PRINTER**



INKJET PRINTER

X-RAY FILMS

KONICA DRY FILM



AGFA FILM




FUJI DIHL DRY FILM



CARESTREAM FILM

SPARE PARTS



Hand Switch

X-RAY TUBE



Bucky

LBD





Excellence Comes Standard



SPOX HEALTHCARE PVT. LTD.

📍 B-15 Sector-6, Noida, Uttar Pradesh 201301 (INDIA)

📞 +91 98113 14143 🌐 www.spox.in

✉ spoxhealthcare@gmail.com

Certified & Approved by:

